
# Quickly, Accurately

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Left: SHRIMP (Sensitive High-Resolution Ion Micro Probe). Upper Right: Laser-ablation ICP-MS. Lower Right: Mineral liberation analyzer.

## Analysis of Critical Metal Ores

To extract critical metals from ores we need to know what minerals and elements are present in the ores. Formation age and isotopic composition are also important in determining the origin of the ores. The Geological Survey of Japan possesses the latest equipment and techniques for quick and accurate analysis of ores, namely laser ablation inductively coupled plasma mass spectrometry (LA-ICPMS) and sensitive high-resolution ion micro-probe (SHRIMP), both of which can determine chemical concentrations and isotopic ratios in small areas with diameters of less than several tens of microns. These methods can be used to detect very small amounts of a wide range of elements, from lithium to uranium, with a sensitivity of less than 1 ppm (0.0001%).

### Highly precise isotope analysis in a micro-area

LA-ICPMS, by instantly vaporizing a small amount of sample and atomizing the elements in a plasma at about 7000 °C, determines chemical composition by mass spectrometry. Analysis of one sample is completed in only in 2 minutes. In contrast, SHRIMP measures the number of secondary ions sputtered from the surface by emitting a primary ion beam to the micro-scale surface area (about 20 μm) under a high vacuum. This instrument is suited particularly to highly precise analyses such as uranium-lead dating and oxygen isotope ratio determination.

The Geological Survey of Japan, by making the best use of these latest technologies and items of equipment, promotes and supports critical metal resource exploration and assessment around the world.

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